

Dynacare Laboratories Test Changes

New

GM1 ANTIBODY

Mnemonic GM1 AB
CPT 83520 x2
Test Code 3550220
Includes
Synonyms
Specimen Serum
Volume 1.0 mL
Container SST or Red Top Tube
Special Inst
Specimen Prep Centrifuge* For Red Top only, transfer serum to Plastic Transport Tube.
Storage Room Temperature
Transport Temp Room Temperature
Method Enzyme-Linked Immunosorbent Assay (ELISA)
Rej Criteria
Reference Value <1:100 titer
Performed One time per week
Comments

Old

GM1 ANTIBODY

Mnemonic GM1 AB
CPT 83520 x2
Test Code 3550220
Includes
Synonyms
Specimen Serum
Volume 1.0 mL
Container SST or Red Top Tube
Special Inst
Specimen Prep Centrifuge* For Red Top only, transfer serum to Plastic Transport Tube.
Storage Room Temperature
Transport Temp Room Temperature
Method Enzyme-Linked Immunosorbent Assay (ELISA)
Rej Criteria
Reference Value Negative;
Positives will be titered
Performed One time per week
Comments

NOTE:

Reference Value Changed

Dynacare Laboratories Test Changes

New

HEPATITIS C RIBA

Mnemonic	HEPC RIBA
CPT	86804
Test Code	3525580
Includes	
Synonyms	
Specimen	Serum
Volume	1.0 mL
Container	SST or Red Top Tube
Special Inst	
Specimen Prep	Centrifuge* For Red Top only, transfer serum to Plastic Transport Tube
Storage	Refrigerate
Transport Temp	Refrigerated
Method	Recombinant Immunoblot Assay (RIBA)
Rej Criteria	
Reference Value	Negative
Performed	Monday - Friday
Comments	

Old

HEPATITIS C RIBA

Mnemonic	HEPC RIBA
CPT	
Test Code	
Includes	
Synonyms	
Specimen	
Volume	
Container	
Special Inst	
Specimen Prep	
Storage	
Transport Temp	Refrigerated
Method	
Rej Criteria	
Reference Value	
Performed	
Comments	Testing is not available due to a national reagent shortage. Clients may consider confirmation testing using Hepatitis C Virus NAAT during this period. Note: this test requires a dedicated sample. This test cannot be added on using the HCV Ab specimen and may require that another specimen be collected on a new test submission

NOTE:

CPT Changed
PDM Changed
Specimen Changed
Volume Changed
Container Changed
Specimen Prep Changed
Storage Changed
Method Changed
Reference Value Changed
Comments Changed

Dynacare Laboratories Test Changes

New

HUMAN HERPES VIRUS TYPE 6 ANTIBODY, IgG

Mnemonic	HHV6 IGG
CPT	86790
Test Code	3550200
Includes	
Synonyms	HHV6, Roseola
Specimen	Serum
Volume	1.0 mL
Container	SST or Red Top Tube
Special Inst	
Specimen Prep	Centrifuge* For Red Top only, transfer serum to Plastic Transport Tube.
Storage	Room Temperature
Transport Temp	Room Temperature
Method	Enzyme immunoassay (EIA)
Rej Criteria	
Reference Value	Negative <0.76 index; Equivocal 0.76 - 0.99 index; Positive >0.99 index
Performed	Two times per week
Comments	

NOTE:

Reference Value Changed

Old

HUMAN HERPES VIRUS TYPE 6 ANTIBODY, IgG

Mnemonic	HHV6 IGG
CPT	86790
Test Code	3550200
Includes	
Synonyms	HHV6, Roseola
Specimen	Serum
Volume	1.0 mL
Container	SST or Red Top Tube
Special Inst	
Specimen Prep	Centrifuge* For Red Top only, transfer serum to Plastic Transport Tube.
Storage	Room Temperature
Transport Temp	Room Temperature
Method	Enzyme immunoassay (EIA)
Rej Criteria	
Reference Value	Negative <0.90 index; Equivocal 0.90 - 1.10 index; Positive >1.10 index
Performed	Two times per week
Comments	

Dynacare Laboratories Test Changes

New

LEUKOCYTE ALKALINE PHOSPHATASE STAIN

Mnemonic	LAP
CPT	85540
Test Code	2001110
Includes	
Synonyms	LAP Score
Specimen	Whole Blood Protected From Light
Volume	7.0 mL
Container	Dark Green Top (Sodium Heparin) Tube
Special Inst	Specimen must be less than 24 hours old.
Specimen Prep	Do NOT centrifuge. 6 unstained slides will be made at Dynacare Laboratory.
Storage	Room Temperature
Transport Temp	Room Temperature
Method	Qualitative Enzymatic
Rej Criteria	
Reference Value	25 - 130
Performed	Monday - Friday
Comments	

NOTE:

Special Instructions Changed
Specimen Prep Changed
Storage Changed
Transport Temp Changed
Reference Value Changed

Old

LEUKOCYTE ALKALINE PHOSPHATASE STAIN

Mnemonic	LAP
CPT	85540
Test Code	2001110
Includes	
Synonyms	LAP Score
Specimen	Whole Blood Protected From Light
Volume	7.0 mL
Container	Dark Green Top (Sodium Heparin) Tube
Special Inst	Specimen must be less than 24 hours old. PROTECT FROM LIGHT
Specimen Prep	Do NOT centrifuge. Transfer to an amber plastic transport tube. If amber tubes are unavailable, cover standard transport tube completely, top and bottom, with aluminum foil. Identify specimen with patient name directly on the container and on the outside of the aluminum foil. Secure with tape.
Storage	Refrigerate and Protect From Light
Transport Temp	Room Temperature. Protect From Light.
Method	Qualitative Enzymatic
Rej Criteria	
Reference Value	Normal: 13 - 130; 0 - 13: Seen in chronic granulocytic leukemia and paroxysmal nocturnal hemoglobinuria. Elevated values with polycythemia vera, myeloid metaplasia and reactive leukocytosis.
Performed	Monday - Friday
Comments	